Technical Data

DU 8009 X-Ray Tube

453570041700

453575083768 Rev A

X-ray Tube Technical Data
5.7 MJ/8.0 MHU CT Insert
0.7 x 0.6/0.9 x 0.9 Focal Spots

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This X-ray Tube will be used as part of a registered and controlled X-Ray Tube Housing Assembly. Please refer to corresponding X-Ray Tube Housing Assembly document for important safety warnings, servicing and operating instructions, warranty and return procedures associated with installation, use and maintenance of this product.

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Copyright

© 2004 DUNLEE

All rights reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright-holder. DUNLEE reserves the right to make changes in specifications or to discontinue any product, at any time without notice or obligation, and is not liable for any consequences resulting from the use of this publication.

Equipment is subject to change without notice. All changes will be in compliance with regulations governing manufacture of medical equipment. Printed in the USA. Document originally drafted in English.

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DISPOSAL

Dunlee Manufactures state-of-the-art x-ray equipment in terms of safety and environmental protection. Under normal operation, there are no risks to persons or the environment. To comply with regulations, it is necessary to use materials that may be harmful to the environment and therefore must be disposed of in a proper manner. Contact your Dunlee sales and service organization for proper handling and disposal.

This metal x-ray tube may contain beryllium. Disposal of this product must be through appropriate channels. (This statement is valid only for metal x-ray tubes, please disregard for glass tubes).

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RECYCLING STATEMENT

Dunlee strives to control manufacturing processes and environmental pollution. Consequently, Dunlee recycles select materials and components from X-ray Tubes. Elaborate testing and quality assurance procedures, including detailed inspection of all subcomponents, ensures the same high level of quality and functionality provided by new materials is maintained.

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INSTALLATION AND ASSEMBLY

The X-ray Tube must be installed into an appropriate x-ray tube housing assembly for safe and proper use. Only trained and qualified personnel must perform this final assembly and be performed in accordance with the documentation and guidelines provided by Dunlee. Failure to follow the appropriate procedure for final assembly will void the warranty.

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

0.04 mm AI @ 70 kV

Applicable Standards
Normes applicables
Anwendbare Standards
Normativa aplicable

IEC 60613, IEC 60552
IEC 60336, EN 60601-1-3
IEC 60878

Section 3:

EMISSION CHARACTERISTICS, CARACTÉRISTIQUES RELATIVES AUX ÉMISSIONS, EMISSIONSMERKMALE, CARACTERÍSTICAS DE EMISIÓN

Section 4:

RATING CHARTS, DIAGRAMMES DE CHARGE, BELASTUNGSKURVEN, CURVAS DE CARGA

MAX. kV AND mA LIMITS

kV	mA 	VCT Hi Power Option	W/O VCT Hi Power
		mA 	mA 
80	300	675	675
100	310	770	700
120	335	800	700
140	335	715	600

453575083769 Rev A

Section 4:

RATING CHARTS, DIAGRAMMES DE CHARGE, BELASTUNGSKURVEN, CURVAS DE CARGA

Single Exposure Limits (mA)

▣ 0.7 x 0.6 ===

Scan Time Sec.	140kV	120kV	100kV	80kV
5	325	335	310	300
10	320	335	310	300
20	305	335	310	300
30	295	335	310	300
40	290	335	310	300
50	280	335	310	300
60	275	335	310	300

453575083771 Rev A

Single Exposure Limits (mA)

■ 0.9 x 0.9 ===


Scan Time Sec.	140kV	120kV	100kV	80kV
5	615*	725*	760*	675
10	575*	680*	760*	675
20	535	630	760	675
30	490	575	695*	675
40	445	525	635	675
50	410	485	585	675
60	385	455	550	675

453575083770 Rev A



*Available with VCT Hi Power Option

Section 5:

-ANODE HEAT AND HEAT DISSIPATION -LA CHALEUR DE L'ANODE ET DISSIPATION THERMIQUE
 -ANODE HITZE UND WÄRMEABLEITUNG -CALOR DEL ÁNODO Y DISIPACIÓN DE CALOR


Section 6: MODULATION TRANSFER FUNCTION, FONCTION DE TRANSFERT DE MODULATION, MODULATIONSÜBERTRAGUNGSFUNKTION, FUNCIÓN DE TRANSFERENCIA DE MODULACIÓN

Section 7:

GENERAL DIMENSIONAL DATA, DONNÉES GÉNÉRALES RELATIVES AUX DIMENSIONS,
ALLGEMEINE MAßANGABEN, DATOS GENERALES SOBRE LAS DIMENSIONES

453575083767 Rev A

Section 8:

Sales, Service, and Manufacturing Locations

Dunlee World Headquarters

555 North Commerce Street
Aurora, Illinois 60504 USA
Tel: 800.238.3780 (U.S. & Canada)
630.585.2100
Fax: 630.585.2125

Dunlee Tubemaster Facility

2312 Avenue J
Arlington, Texas 76006 USA
Tel: 800.544.9729 (U.S. & Canada)
817.640.7666
Fax: 817.640.6644

Dunlee Medical Components

Sales and Service Office - Germany
Marbacher Strasse 114
D-40597 Düsseldorf, Germany
Tel: 49.0.211.20.30.05
Fax: 49.0.211.20.30.00

Dunlee Richmond Facility

8819 Whitepine Road
Richmond, Virginia 23237 USA
Tel: 800.526.0555 (U.S. & Canada)
804.714.2501
Fax: 804.714.2507

Philips Sales Office - Mexico

La Palma 6, 3rd Floor
Col. San Fernando La Herradura
Huixquilucan, Edo. Mex., México.
CP 52736
Tel: 52.55.52.69.92.03
52.55.21.09.75.21
Fax: 52.55.52.69.92.04

Philips Sales Office - Europe

Philips Medical Systems S.p.A
Via Casati 23
20052 Monza (MI) Italy
Tel: 39.039.203.6536
Fax: 39.039.203.6175

Philips Sales Office - Asia

Level 9, Three Pacific Place
1 Queen's Road East
Wan Chai
Hong Kong
Tel: 630.585.2100
Fax: 630.585.2505

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European Union Representative

Philips Medical Systems
Dept.: Corporate Quality & Regulatory
Veenpluis 4
5684 PC Best
The Netherlands

DUNLEE

Visit us at www.dunlee.com

DUNLEE