
DUNLEE

INSTRUCTIONS FOR USE

X-RAY TUBE ASSEMBLIES
3RDPB

B Ξ YOND SOLUTIONS

Table of Contents

English, English (US).....	5
Deutsch, German.....	25
Français, French.....	45
Español, Spanish.....	67
体中文, Chinese simplified.....	89
العربية, Arabic.....	107
Български, Bulgarian.....	125
Hrvatski, Croatian.....	147
Česky, Czech.....	167
Dansk, Danish.....	187
Nederlands, Dutch.....	207
Eesti, Estonian.....	229
Suomi, Finnish.....	249
Ελληνικά, Greek.....	269
Magyar, Hungarian.....	291
Íslenska, Icelandic.....	311
Italiano, Italian.....	331
日本語, Japanese.....	351
한국어, Korean.....	371
Latviešu valodā, Latvian.....	391
Lietuvių, Lithuanian.....	411

Table of Contents

Malti, Maltese.....	431
Norsk, Norwegian.....	451
Polski, Polish.....	471
Português, Portuguese.....	491
Română, Romanian.....	511
Русский, Russian.....	531
Srpski, Serbian.....	553
Slovensky, Slovak.....	573
Slovenščina, Slovenian.....	593
Svenska, Swedish.....	613
Türkçe, Turkish.....	633
Українська, Ukrainian.....	653

Instructions for Use

X-Ray Tube Assemblies

English, English (US)

Table of contents

1 Document Information.....	6
2 Intended Purpose.....	7
3 Safety Information.....	8
Safety Message Used in this Document.....	8
General Safety Information.....	8
Limits for Operation, Storage, and Transportation.....	9
Electrical Safety.....	11
Radiation Protection.....	11
Thermal Safety.....	11
Cooling / Insulating Media.....	12
Electromagnetic Compatibility (EMC).....	12
Decommissioning.....	13
China RoHS Label and Material Declaration Table (China only).....	15
Ukraine RoHS Label and Representative (Ukraine only).....	16
4 Conformity.....	17
5 Compatibility.....	18
6 Maintenance.....	20
Planned Maintenance by the Organization that is Responsible for the Operation.....	20
Cleaning and Disinfection.....	21
Planned Maintenance Schedule.....	21
7 Product Identification Label.....	23

1 Document Information

Document data

Doc-ID: 300004427491C

Release date: 03/2020

Copyright

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Disclaimer

All rights reserved. Koninklijke Philips N.V. reserves the right to make changes in specifications in accordance with the terms in commercial or other applicable agreements between supplier and buyer, and will not be liable for any consequences resulting from the use of this publication. Reproduction in whole or in part is prohibited without the prior written consent of the copyright owner. The English version of this document is the original version.

Product tradename

The duplicate product identification label which accompanies this Instructions For Use provides the tradename of the X-ray tube assembly.

Contact information of the manufacturer

Philips Medical Systems DMC GmbH

Röntgenstraße 24

22335 Hamburg / Germany

For further help see different contact options on: www.dunlee.com

2 Intended Purpose

Intended use

The X-ray tube assembly is an X-ray generating tube encased in a radiation-shielded housing. The X-ray tube assembly is intended to emit X-ray photons to enable the creation of diagnostic images of human patients. The X-ray tube assembly is only used as integrated part of an X-ray system and has no medical purpose by its own. The intended treatment, duration and treatment parameters are not defined on X-ray tube assembly level. This is defined by the manufacturer of the X-ray system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the X-ray system.

Intended user population

The X-ray tube assembly is remotely controlled by the X-ray system, which itself is operated by qualified clinical users only. Further details regarding the intended user and patient population are defined by the manufacturer of the X-ray system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the X-ray system. The X-ray tube assembly must be installed, repaired and maintained only by technically qualified field service engineers who have received the appropriate instructions and training and who are authorized in particular for the installation and commissioning of radiography equipment in medical applications.

Intended context

The X-ray tube assembly is intended to be used in a non-condensing, climate controlled, indoor, clinical environment and is not intended to be used for homecare or by lay persons. The X-ray tube assembly is intended to be used in permanently installed stationary X-ray systems. The use frequency is not limited within the expected service lifetime. This X-ray tube assembly is intended to be re-used. The re-use of this X-ray tube assembly does not require any special treatment. Additional treatment for re-use may be required and verified by the manufacturer of the X-ray system for specific applications.

Operating principles

Contraindications are not defined on X-ray tube assembly level, because this is defined by the manufacturer of the X-ray system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the X-ray system. Specific warnings are provided within the accompanying documents of the X-ray tube assembly. The expected service lifetime is referenced in the risk management file and stated in the accompanying documents of the X-ray tube assembly.

Essential performance

The X-ray tube assembly does neither have essential performance per se. Nor is any function of the X-ray tube assembly intended to contribute to the essential performance of the X-ray system in which it is integrated.

3 Safety Information

Warning

These Instructions For Use are made to make it possible to work safely with the X-ray tube assembly. Operate the X-ray tube assembly only in compliance with the instructions in this document and the operating instructions of the applying X-ray system. Do not use it for purposes other than for that it is intended.

Read me,
understand me,
and obey me!

Safety Message Used in this Document

This safety-relevant message occur in this instruction for use:

WARNING

shows a hazardous situation for people. If you do not prevent this situation, there is a risk of death or serious injury.

General Safety Information

The for the operation responsible organization of the applying X-ray system is always responsible for the compliance with the regulations that apply to the installation and the operation of this X-ray tube assembly.

- Only technically qualified field service engineers who have received the appropriate instructions and training for equipment in medical applications and in particular for the installation and commissioning of this X-ray tube assembly are allowed to install and repair this equipment. The installation, calibration and testing of the X-ray tube assembly must be performed in accordance with the respective instructions of the applying X-ray system. Information to verify that the X-ray tube assembly is properly installed and it is ready to perform safely as intended, shall be documented in the IEC 62353 Test Report. The IEC 62353 Test Report is part of the accompanying documents of this X-ray tube assembly. It has to be filled in and to be signed at the end of the installation. The signed IEC 62353 Test Report has to be filed along with the reference documents of the applying X-ray system.

- Only trained and authorized personal are permitted to operate the X-ray tube assembly.
- If the X-ray tube assembly has electrical, mechanical, or functional defects, do not use it. This fact particularly applies to faulty indicators, displays, warnings, and alarms.
- If functional defects or other deviations from the normal operational behavior occur, set the X-ray system to OFF. Tell the company that puts the X-ray system that includes this X-ray tube assembly on the market about the defects or deviations immediately. Continue the operation of the X-ray system only after it is repaired. An operation with defective components increases safety risks.
- Do not install the X-ray tube assembly to other medical devices than this X-ray tube assembly is made for.
- The manufacturer of the X-ray tube assembly is responsible for the safety features of this X-ray tube assembly. The manufacturer rejects the responsibility for the safety features of the X-ray tube assembly if modifications are made.
- For the X-ray tube assembly correct operation and regular, competent maintenance are necessary.
- Operate the X-ray tube assembly correctly. Maintain or let maintain the X-ray tube assembly correctly. In other cases, the manufacturer of the X-ray tube assembly is not responsible for incorrect operation, damages, or injuries that occur.
- Do not make or let make service or maintenance at the X-ray tube assembly while the X-ray tube assembly is in use.
- Any serious incident that occurs in relation to the X-ray tube assembly, must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

300004427491C / D000411798D / 740 * 03/2020

Limits for Operation, Storage, and Transportation

Safe operation of the X-ray tube assembly is only made sure when you use the X-ray tube assembly within the limits of its specifications.

If you ignore the operation limits, there is the danger of coolant that flows out. The X-ray tube assembly becomes too hot. Component parts that are too hot can blow up as a result of an implosion or explosion. Refer to the respective operation limits of the applying X-ray system.

- Do not operate the X-ray tube assembly in presence of a flammable anesthetic mixture with air or oxygen or a nitrous oxide.
- Do not use detergents and disinfectants, also the detergents and disinfectants used on the patient, that can create explosive gas mixtures.
- The oxygen content of the ambient air during operation must be less than 25 %.

- Obey the specified environmental limits for the X-ray tube assembly:


Environmental limits	X-ray tube assembly		
	REEVO240G	DA200ULTRA	DA200P40
In operation			
Ambient temperature [°C]			
Minimum	+15	+5	+18
Maximum	+35	+35	+26
Ambient humidity [%]			
non-condensing			
Minimum	+20	+30	+30
Maximum	+80	+60	+60
Ambient atmospheric pressure [kPa]			
Minimum	+70	+70	+70
about 3,000 m above the sea level			
Maximum	+110	+110	+110
about -750 m below the sea level			
For storage and transportation			
Temperature limits [°C]			
Minimum	-25	-25	-25
Maximum	+70	+70	+70
Humidity limits [%]			
Minimum	---	---	---
Maximum	+95	+95	+95

300004427491C / D000411798D/740 * 03/2020

Dunlee

Electrical Safety

Use this X-ray tube assembly only in medical rooms that obey the requirements of relevant local or transregional standards and laws.

Warning:
Electricity

Radiation Protection

Before you make an X-ray exposure, put all necessary radiation precautions into effect. You can find information about radiation precautions in the instructions for use of the applying X-ray system.

- Before you start with radiation, make sure that the X-ray system supplies sufficient proof of a correct image quality.
- To stop radiation, use the emergency stop switch of the X-ray system!

Warning:
Ionizing radiation

Thermal Safety

To prevent damage due to thermal overloading, make sure that you do not operate the X-ray tube assembly outside its specified load parameters. This prevention decreases the risks to the patient, operating staff, third parties, and the environment.

Warning:
Explosive material

Cooling / Insulating Media

- Do not swallow media that spills.
- Do not let media or its vapors to go into your sewerage system.
- Do not breathe the media vapors.
- Make sure to supply closed rooms with fresh air.
- Do not cause vomiting.
- Remove spilled media. Use a liquid-absorbent material.
- Discard the cleaning materials. Obey the local environmental laws and regulations.

Warning:
General

Electromagnetic Compatibility (EMC)

In accordance with its intended use, this electronic product obeys the law that controls EMC. This law gives the permitted emission levels from electronic equipment and its necessary stability against electromagnetic fields.

Electronic products that obey the EMC requirements have a design so that under normal conditions there is no risk of incorrect operation because of electromagnetic interference. If radio signals from high-frequency transmitters with a high transmitting power occur, the risk of electromagnetic incompatibility when operated closely to medical devices can occur.

In unusual circumstances, unintended functions of the electronic product can occur. Possibly undesirable risks for the patient, operating staff, or third parties occur. For this reason, prevent all kinds of transmission with mobile devices. The risks also apply when the medical device is in STANDBY mode.

The X-ray tube assembly is suitable for professional healthcare facility environment.

Too high EMC disturbances can lead to unwanted X-ray, interrupted X-ray, wrong X-ray, unexpected status changes, error messages, component defects, or loss of installation parameters.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the X-ray tube assembly. Otherwise, degradation of the performance of this equipment could result.

The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The for the operation responsible organization might need to take mitigation measures, such as relocating or re-orienting the equipment.

Warning:
Non-ionizing radia-
tion

Decommissioning

The take-back obligation, correct disposal, and recovery of the X-ray tube assembly refer to the European Waste Electrical and Electronic Equipment (WEEE) directive. It also refers to the requirements of local and transregional legal requirements.

The manufacturer of the X-ray tube assembly assembles state-of-the-art X-ray tube assemblies in terms of safety and environmental protection. If no parts of the X-ray tube assembly are opened and if the X-ray tube assembly is used correctly there are no risks to persons or the environment.

To obey regulations, sometimes it is necessary to use materials that are harmful to the environment. Discard these materials in a correct manner.

This X-ray tube assembly contains materials that are toxic. Do not discard the X-ray tube assembly together with industrial or domestic waste.

Obey the disposal regulations!

The manufacturer

- supports you in the disposal of the X-ray tube assembly in accordance with valid legal requirements.
- takes back the X-ray tube assembly.
- returns re-usable parts to the production cycle via certified disposal companies. Extensive test and quality assurance procedures as well as detailed checks of the components guarantee the same high level concerning quality and functionality that is expected from new materials.
- makes a contribution to the protection of the environment.

In case you have questions concerning safe disposal, please consult the manufacturer in full confidence.

China RoHS Label and Material Declaration Table (China only)

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Ukraine RoHS Label and Representative (Ukraine only)

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Philips Medical Systems DMC GmbH
Röntgenstraße 24, 22335 Hamburg, Germany

Ukraine representative for X-ray tubes, which are marked with the conformity mark of Ukraine RoHS Technical Regulation:

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Dunlee

4 Conformity

The X-ray tube assembly that is shown in these instructions for use puts the provisions of the

- European CE Marking of Conformity, and

Identification num-
ber of the notified
body

- FDA 21 CFR 1020.30

into effect.

If you have further questions that refer to local or transregional legal requirements,
make a contact to:

Philips Medical Systems DMC GmbH
Röntgenstraße 24
22335 Hamburg / Germany

For further help see different contact options on: <http://www.dunlee.com>

5 Compatibility

This X-ray tube assembly is designed and manufactured to be compatible with specific diagnostic imaging equipment. Prior to installation, make sure that the X-ray tube assembly and the applying X-ray system are compatible.

The matrix below identifies compatible X-ray systems with an "X".

X-ray systems which do not have an "X" or are not listed are not validated.

X-ray system	X-ray tube assembly		
	REEVO240G	DA200ULTRA	DA200P40
GE BrightSpeed™ Edge		X	
GE BrightSpeed™ Elite		X	
GE BrightSpeed™ Excel		X	
GE Discovery™ CT590 RT	X		
GE Discovery™ NM/CT570c	X		
GE Discovery™ NM/CT 670		X	
GE Discovery™ NM/CT 670 Pro		X	
GE Discovery™ PET/CT 610		X	
GE Discovery™ PET/CT 710		X	
GE Discovery™ RX		X	
GE Discovery™ ST		X	
GE Discovery™ STE		X	
GE Discovery™ VCT	X		
GE HiSpeed™ QX/i		X	
GE LightSpeed™ 16		X	
GE LightSpeed™ Plus		X	
GE LightSpeed™ Pro 16	X		
GE LightSpeed™ Pro 32	X		
GE LightSpeed™ QXI		X	
GE LightSpeed™ RT16	X		
GE LightSpeed™ Ultra		X	
GE LightSpeed™ VCT	X		
GE LightSpeed™ VCT Select	X		
GE LightSpeed™ VCT XT	X		

X-ray system	X-ray tube assembly		
	REEVO240G	DA200ULTRA	DA200P40
GE LightSpeed™ VCT XTe	X		
GE LightSpeed™ Xtra	X		
GE Optima™ CT540		X	
GE Optima™ CT580	X		
GE Optima™ CT580 RT	X		
GE Optima™ CT580 W	X		
GE Optima™ CT 660			X
GE Optima™ PET/CT 560		X	
GE Optima™ PET/CT 560 FX		X	

6 Maintenance

Only technically qualified field service engineers who have received the appropriate instructions and training for the equipment in medical applications and in particular for the installation and commissioning of this X-ray tube assembly are allowed to install and repair this equipment. The installation, calibration and testing of the X-ray tube assembly must be performed in accordance with the respective instructions of the applying X-ray system.

Information to verify that the X-Ray tube assembly is properly installed and it is ready to perform safely as intended, shall be documented in the IEC 62353 Test Report. The IEC 62353 Test Report is part of the accompanying documents of this X-Ray tube assembly. It has to be filled in at the end of the installation.

The signed IEC 62353 Test Report has to be filed along with the reference documents of the applying X-Ray system.

Planned Maintenance by the Organization that is Responsible for the Operation

The for the operation responsible organization of medical devices is committed to do planned maintenance. This planned maintenance must refer to the regulations for the prevention of industrial accidents, local or transregional medical product laws, and further regulations.

This X-ray tube assembly requires regular planned maintenance:

Make sure that service personnel do planned maintenance on this X-ray tube assembly regularly.

Precautionary planned maintenance measures prevent personal injury and make sure that the for the operation responsible organization puts all commitments into effect.

If functional defects or other deviations from the normal operational behavior occur, set the X-ray system to OFF. Tell the company that puts the X-ray system that includes this X-ray tube assembly on the market about the defects or deviations immediately. Continue the operation of the X-ray system only after it is repaired. An operation with defective components increases safety risks.

Regular maintenance cycles decrease the probability of a spontaneous failure. You cannot prevent spontaneous failures completely. If it is necessary to have the guaranteed availability of the X-ray system, make sure that a redundant standby X-ray system is available.

The structural integrity of the X-ray tube assembly does not wear.

Cleaning and Disinfection

Hazard of overdosed radiation or electrical shock.

- Set the X-ray system to OFF before you start the cleaning procedure of the X-ray tube assembly.
- Do not open the covers of the X-ray system.
- Do not change or remove the covers or other parts of the X-ray tube assembly.
The X-ray tube assembly must be directly accessible.
- Do not change the configuration of the X-ray tube assembly.
- Make sure that no detergents get inside the X-ray equipment.
- If you use cleaning agents such as detergents, make sure that they do not contain explosive substances as they can create explosive gas mixtures.
- Obey the regional legal requirements that refer to the cleaning procedure.
- Obey the manufacturers information and safety instructions of the detergents.

If you do not obey this instruction, there is a risk of death or serious injury.

Cleaning

Do not clean the X-ray tube assembly. It is behind the covers of the X-ray system.
A cleaning of the X-ray tube assembly is not necessary.

Disinfection

A disinfection of the X-ray tube assembly is not necessary, unless the system manufacturer requires it.

Planned Maintenance Schedule

- Do checks of the X-ray tube assembly for apparent defects. Do the checks unless they are contrary to the configuration of the X-ray system.


Interval	Scope of work
Always you operate the X-ray system	<p>Do a check if error messages of the X-ray system occur. If the X-ray tube assembly is thermally overloaded and does not operate, wait until a cooling period lets you operate the X-ray tube assembly again.</p> <p>Do a check for visible media leakage or other contamination.</p>
Weekly	Do a check for unusual noises.

Interval	Scope of work
As per relevant local or transregional standards and laws	Do a constancy test.
As per relevant local or transregional standards and laws	Do a check for image quality.

- If defects occur that refer to the X-ray tube assembly, make a contact to the company that puts the X-ray system that includes this X-ray tube assembly on the market.

7 Product Identification Label

The performance characteristics of the device are defined on the product identification label as part of the accompanying documents.

No.	Label	Designation
01		Registered trade name / trade mark of the manufacturer
02	Philips Medical Systems DMC GmbH Röntgenstraße 24 22335 Hamburg / GERMANY	Address of the legal manufacturer
03	 yyyy-mm	Address of the legal manufacturer yyyy: year of manufacture mm: month of manufacture
04	X-RAY TUBE HOUSING ASSEMBLY	Type of the medical device: X-ray tube housing assembly
05	Model	Name of the product
06	REF	Type number of the product
07	SN	Serial Number of the product
08	MD	Medical device product
09		Value of the permanent filtration
10		Dimensions of the small focus and the respective standard
11		Dimensions of the large focus and the respective standard
12	Nominal X-ray Tube Voltage	Value of the nominal X-ray tube voltage and the respective standard
13	X-RAY TUBE	X-ray tube
14	US only	Compliance statement, for US only
THIS PRODUCT COMPLIES WITH THE DHHS REQUIREMENT OF 21 CFR SUB-CHAPTER J		

No.	Label	Designation
15	MANUFACTURED: month year	Month and year of manufacture
16		The Unique Device Identification [UDI] Data Matrix Code consists of the Global Trade Item Number [GTIN], Serial Number [SN], and Type Number [TN] of the product.
17	(01)xxxxxxxxxxxx	Global Trade Item Number [GTIN]
18	(21)xxxxxxxxxxxx	Serial Number [SN] of the product xxxxxwyyyyy xxxxx : X-ray tube housing w : X-ray tube housing type yyyyy : X-ray tube
19		CSA symbol
20		CE symbol and the number of the notified body
21		Disposal symbol
22		General warning symbol
23		Electronic version of the Instructions For Use. See: www.dunlee.com/IFU
24		China RoHS symbol
25	Duplicate label	Duplicate label for the label overview of the X-ray system
26	(CPO) Refurbished	Certified pre-owned refurbished product
27		Additional label: Instruction manual/booklet must be read
28	 UA.TR.116	Ukraine RoHS Technical Regulation symbol

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All rights reserved. Koninklijke Philips N.V. は、納入業者と購入者の間の商業契約またはその他の適用される契約の条件に従って仕様を変更する権利を留保します。また、この印刷物を使用したことの結果について一切の責任を負いません。著作権保有者から事前に書面による許可を得ることなく本書を複製・転写することは、その全部或いは一部に限らず禁じます。本書の英語版が原本です。