

# CT CROSS REFERENCE

07/01/2018

Subject to change without notice. All products listed may be trademarked by the associated OEM.

Brand	CT Scanner Model	Dunlee Replacement Tube	Dunlee Part #
GE	BrightSpeed™ Edge	DA200 Ultra	9896-055-88704
GE	BrightSpeed™ Elite	DA200 Ultra	9896-055-88704
GE	BrightSpeed™ Excel	DA200 Ultra	9896-055-88704
GE	BrightSpeed Edge™ (Select / Lite)	DA165NP	9896-055-83901
GE	BrightSpeed Elite™ (Select / Lite)	DA165NP	9896-055-83901
GE	BrightSpeed Excel™ (Select / Lite)	DA165NP	9896-055-83901
GE	Brivo™ CT315/325	DA135Ct/e	9896-055-80800
GE	Brivo™ CT385	DA135Ct/e	9896-055-80800
GE	Discovery NM/CT570c	Reevo 240 G	9896-055-96401
GE	Discovery™ CT590 RT	Reevo 240 G	9896-055-96401
GE	Discovery™ NM/CT 670	DA200 Ultra	9896-055-88704
GE	Discovery™ NM/CT 670 Pro	DA200 Ultra	9896-055-88704
GE	Discovery™ PET/CT 610	DA200 Ultra	9896-055-88704
GE	Discovery™ PET/CT 710	DA200 Ultra	9896-055-88704
GE	Discovery™ RX	DA200 Ultra	9896-055-88704
GE	Discovery™ ST (PET/CT)	DA200 Ultra	9896-055-88704
GE	Discovery™ STE (PET/CT)	DA200 Ultra	9896-055-88704
GE	Discovery™ VCT	Reevo 240 G	9896-055-96401
GE	HiSpeed CT/e, CT Dual™	DA135Ct/e	9896-055-80800
GE	HiSpeed DX/i	DA165NP	9896-055-83901
GE	HiSpeed FX/i	DA165NP	9896-055-83901
GE	HiSpeed LX/i	DA165NP	9896-055-83901



Brand	CT Scanner Model	Dunlee Replacement Tube	Dunlee Part #
GE	LightSpeed™ 16	DA200 Ultra	9896-055-88704
GE	LightSpeed™ Plus	DA200 Ultra	9896-055-88704
GE	LightSpeed™ Pro <sup>32</sup>	Reevo 240 G	9896-055-96401
GE	LightSpeed™ QX/i	DA200 Ultra	9896-055-88704
GE	LightSpeed™ RT <sup>16</sup>	Reevo 240 G	9896-055-96401
GE	LightSpeed™ Ultra	DA200 Ultra	9896-055-88704
GE	LightSpeed™ VCT	Reevo 240 G	9896-055-96401
GE	LightSpeed™ VCT Select	Reevo 240 G	9896-055-96401
GE	LightSpeed™ VCT XT	Reevo 240 G	9896-055-96401
GE	LightSpeed™ VCT XTe	Reevo 240 G	9896-055-96401
GE	LightSpeed™ Xtra	Reevo 240 G	9896-055-96401
GE	Optima™ CT520	DA165NP	9896-055-83901
GE	Optima™ CT520 Pro	DA165NP	9896-055-83901
GE	Optima™ CT540	DA200 Ultra	9896-055-88704
GE	Optima™ CT580	Reevo 240 G	9896-055-96401
GE	Optima™ CT580 RT	Reevo 240 G	9896-055-96401
GE	Optima™ CT580 W	Reevo 240 G	9896-055-96401
GE	Optima™ CT660	DA200 P40	4598-006-53451
GE	Optima™ PET/CT 560	DA200 Ultra	9896-055-88704
GE	Optima™ PET/CT 560 FX	DA200 Ultra	9896-055-88704
Siemens	Somatom Sensation 4	S 532B/Akron B	9896-055-95702
Siemens	Somatom Sensation Volume Zoom	S 532B/Akron B	9896-055-95702
Siemens	Somatom Sensation Volume Access	S 532B/Akron B	9896-055-95702
Siemens	Somatom Sensation 10	S 532Q/Akron Q	9896-055-95902
Siemens	Somatom Sensation 16	S 532Q/Akron Q	9896-055-95902
Siemens	Somatom Sensation Cardiac	S 532Q/Akron Q	9896-055-95902

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Dunlee reserves the right to make changes in specification and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Dunlee is a brand of the Philips Company Group

**Contact us for further information and to find the tube best suited to your needs.**

Philips Medical Systems DMC GmbH  
Tel.: +49 40 5078-1391 • Fax: +49 40 5078-1611  
Röntgenstrasse 24 • 22335 Hamburg • Germany

**dunlee.com**